



PRENATAL REFERENCE SHEET

TUBE TOP COLOR	# OF TUBES	TESTS MENU	SEND TO	TURN AROUND TIME FROM RECEIPT OF SAMPLE
	1 per	1052 <input type="checkbox"/> Fetal Fibronectin	LENETIX [®]	24 - 72 Hours • Vaginal Swab of patient • Refrigerate until shipment
	1 per	<input type="checkbox"/> Paternity Testing	LENETIX [®]	3 - 5 Days • Swab firmly inside of each patient's mouth area. One swab kit per person. • Refrigerate until shipment
	1	<p>FIRST & SECOND TRIMESTER <input type="checkbox"/> Integrated 190 Test 3327 <input type="checkbox"/> Stage 1: PAPP-A (11 weeks 0 days - 13 weeks 6 days) 3325 <input type="checkbox"/> Stage 2: AFP, total βhCG, uE3, inhibin A (15 weeks 0 days - 22 weeks 6 days)</p> <p><input type="checkbox"/> Modified Sequential 30•190 <input type="checkbox"/> Contingent 30•1500•190 1414 <input type="checkbox"/> Stage 1: PAPP-A, total βhCG (11 weeks 0 days - 13 weeks 6 days) 3325 <input type="checkbox"/> Stage 2: AFP, total βhCG, uE3, inhibin A (15 weeks 0 days - 22 weeks 6 days)</p> <p>FIRST TRIMESTER - (11 weeks 0 days - 13 weeks 6 days) 1414 <input type="checkbox"/> Combined (total βhCG, PAPP-A) 1414A <input type="checkbox"/> Combined Plus (total βhCG, PAPP-A, inhibin A)</p> <p>SECOND TRIMESTER - (15 weeks 0 days - 22 weeks 6 days) 3121 <input type="checkbox"/> Quad Screen (AFP, total βhCG, uE3, inhibin A) 1520 <input type="checkbox"/> Triple Screen (AFP, total βhCG, uE3) 608 <input type="checkbox"/> AFP (ONTD Only) 1101 <input type="checkbox"/> Inhibin A 0505 <input type="checkbox"/> Inhibin B <input type="checkbox"/> Other</p>	LENETIX [®]	<p>Stage 1 Held until Stage 2</p> <p>Stage 2 24 - 72 hours</p> <p>All others 24 - 72 hours</p> <p>Handling: • Allow to clot about 30 minutes and spin 2500 rpm for 10 minutes. • Refrigerate until shipment</p>
	1-3	<p>1655, 3108 <input type="checkbox"/> Amniotic Fluid Chromosome Analysis (Includes Amniotic Fluid AFP) 1118 <input type="checkbox"/> Chorionic Villus Sampling (CVS) 1116 <input type="checkbox"/> POC/Other Tissue, Specify: _____</p> <hr/> <p>1094 <input type="checkbox"/> AneuVysion[®] FISH (Chromosomes 13, 18, 21, X, Y) 3108 <input type="checkbox"/> Amniotic Fluid AFP</p> <hr/> <p>82103 <input type="checkbox"/> Acetylcholinesterase (ACHE)</p>	LENETIX [®] Handling: • Maintain room temperature • DO NOT FREEZE	<p>7 - 15 days</p> <hr/> <p>24 - 48 hours</p> <hr/> <p>3 to 5 days</p>
	1	<p>684 <input type="checkbox"/> Peripheral Blood Chromosome Analysis</p> <hr/> <p>1094 <input type="checkbox"/> AneuVysion[®] FISH (Chromosomes 13, 18, 21, X, Y)</p>	LENETIX [®]	<p>Handling: 10 - 30 days • See box below</p> <hr/> <p>Handling: 24 - 48 hours • See box below</p>
	4	556 <input type="checkbox"/> Tay-Sachs (Enzyme)	LENETIX [®]	<p>Handling: 7 - 10 days • Draw 4 full tubes • Invert tubes 8 times to prevent clotting • Refrigerate until shipment</p>
	2-5	<p>3417 <input type="checkbox"/> Ashkenazi Jewish Carrier 9 Risk Assessment: Cystic Fibrosis, Tay-Sachs, Canavan Disease, Gaucher Disease, Bloom Syndrome, Niemann-Pick Disease (Type A and Type B), Fanconi Anemia (Group C), Familial Dysautonomia, Mucopolidosis (Type IV)</p> <p>1278 <input type="checkbox"/> Bloom Syndrome 1508 <input type="checkbox"/> Canavan Disease 1800 <input type="checkbox"/> CF (5T-Allele) 1341 <input type="checkbox"/> Cystic Fibrosis 3135 <input type="checkbox"/> Familial Dysautonomia 1788 <input type="checkbox"/> Familial Hyperinsulinism</p> <p>738 <input type="checkbox"/> Fanconi Anemia (Group C) 2006 <input type="checkbox"/> Fragile X Syndrome 1509 <input type="checkbox"/> Gaucher Disease 318 <input type="checkbox"/> Glycogen Storage (Type 1A) 319 <input type="checkbox"/> Maple Syrup Urine Disease 6410 <input type="checkbox"/> Mucopolidosis (Type IV) 3127 <input type="checkbox"/> Niemann-Pick Disease (Type A & B) 243 <input type="checkbox"/> Sickle Cell Anemia 558 <input type="checkbox"/> Tay-Sachs DNA (Reflex) 2313 <input type="checkbox"/> Y Chromosome Microdeletion <input type="checkbox"/> Other _____</p>	LENETIX [®]	<p>Call LENETIX[®] for specific turn around times</p> <p>Handling: • Draw 2 - 5 full tubes depending on number of tests ordered. • Invert tubes 8 times to prevent clotting • Refrigerate until shipment</p>
	3	2312 <input type="checkbox"/> SRY, X & Y	LENETIX [®]	<p>Handling: 3 - 5 Days • Draw 3 White and 1 Purple tubes. • Invert 8 times to prevent clotting • Refrigerate until shipment</p>
	1	1118 <input type="checkbox"/> RhD & SRY Genotyping (must be drawn after 15 weeks gestational age)	LENETIX [®]	<p>MUST DRAW & SHIP MONDAY-THURSDAY ONLY! • DO NOT FREEZE</p>